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RECOSYN®

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RECOSYN®

1. Intended Use

The product is a medical device used in humans for intra-articular injection in the knee. The device is intended to reduce the functional impairment, to ease the pain and thus to improve quality of life.

2. Indication

It is indicated for patients with symptomatic femorotibial knee osteoarthritis.

3. Description

The product consists of a sterile, isotonic, viscoelastic solution for intra-articular use, supplied in a ready-to-use syringe. One ready-to-use syringe contains 2 ml of the viscoelastic solution. 1 ml of the product contains 10 mg (1 %) of sodium hyaluronate together with sodium chloride, sodium monohydrogen phosphate, citric acid and water for injection. The product is provided sterile and is intended for single use only. The product is sterilized using steam within the blister. No materials of animal origin are used during manufacturing of the product or as raw materials. No pharmaceuticals are integrated into the product.

4. Method of Administration and Dosage

Take the ready-to-use syringe out of the blister, remove the rubber cap from the Luer-Lock adapter by turning it gently. A 21G needle is recommended to be used for intra-articular injections. Attach the needle to the Luer-Lock adapter and lock it into place by turning carefully. During administration, the syringe should be held as shown in figure 1. It is recommended to administer three injections of 2 ml of the product per knee at weekly intervals. Both knees may be treated simultaneously and treatment cycles up to three injections may be repeated. In order to avoid an intra-articular infection, a strict aseptic injection technique must be applied. After administration of the product, it is recommended that an ice-pack be placed on the treated joint for 5–10 minutes in order to prevent pain, swelling and bleeding.

5. Performance and Mode of Action

In patients with degenerative joint disease (osteoarthritis), the viscoelasticity of the synovial fluid is significantly impaired. Mechanical stress on the joint leads to the breakdown of the articular cartilage resulting in restricted and painful joint movement. The lubricating properties of this product help reduce pain and improve joint mobility. The performance of the device lasts up to six months following a treatment cycle of three intra-articular injections.

6. Contraindications and Possible Interactions

The device must not be used in patients: with hemochromatosis, ochronosis or hemophilia I with Paget's disease of bone, chondromatosis or villonodular synovitis I with inflammatory, infectious or metabolic arthritis (rheumatoid arthritis, ankylosing spondylitis, abarticular rheumatism...) I with history of homeostasis disorders or under anticoagulant treatment I with known allergy to the product or one of its constituents I requiring regular or intermittent steroid therapy I immunocompromised patient or patient presenting with a serious or progressive disease (cardiac, pulmonary, hepatic, renal, hematological, neoplastic or infectious disease) I with venous or lymphatic stasis of the target limb I with skin diseases or infections in the area of the injection site I with a severe acute or chronic disease which the doctor deems incompatible with the product or if the application of the product exposes the patients to additional risk

Currently there is no information of adverse interactions with other intra-articular treatments.

7. Adverse Reactions and Adverse Events

Doctors must inform the patients that there are potential side effects and/or incompatibilities associated with implantation of this device, which may occur immediately or may be delayed. After administration, patients may experience local symptoms in the joint being treated (pain, sensation of heat, reddening and swelling). The following adverse events have been reported for similar products: mild or moderate arthralgia, in rare cases skin rash, pseudoseptic joint effusion, pruritus and muscular cramps. Further adverse events that have been observed in very rare cases are: allergic reactions, anaphylactic shock, hemarthrosis, phlebitis, severe acute inflammatory reaction (SAIR), nasopharyngitis, joint stiffness, tendonitis, bursitis, fever, myalgia and septic arthritis. Doctors must inform the patients to consult their doctor as soon as possible in case of appearance of the adverse events. The doctor must treat these side effects appropriately. Any undesirable effects associated with injection of the product must be reported to the distributor and/or the manufacturer.

8. Warnings

The product must be administered only by a doctor trained to apply intra-articular injections. He or she should also be familiar with all the immunological and other potential risks associated with the use of biological material. The product must not be used in pregnant or breastfeeding women and in patients under 18 years of age since no clinical data is available. Keep the product out of sight and reach of children. The product is intended for single use and the syringe must not in any case be re-sterilized. The re-use of the product creates a potential infection risk for patients or users. Do not use if the blister packaging has been opened and/or is damaged. Do not use a syringe with an open or shifted tip cap in the blister. Do not administer after the expiry date. Patients should be advised to relative rest (but no immobilisation) for 24 hours after each injection in order to avoid strain on treated joints. Discard the syringe, remaining product and needle after use in a special container. A 21G needle is recommended to inject the product. The compatibility of other needles/cannulas with the product and the syringe system has not been tested.

9. Storage

The product should be stored at room temperature (2–25 °C/36–77 °F) and protected from frost and moisture.

Available package sizes: 1, 5 ready-to-use syringes

10. Last revised

2020–01

Version Code: Z07

Manufacturer

CROMA-PHARMA GmbH

Industriezeile 6

2100 Leobendorf

Austria

Tel.: +43 2262 684 68–0

Fax: +43 2262 684 68–15

Distributor (in the Federal Republic of Germany):

Recordati Pharma GmbH

Eberhard-Finckh-Straße 55


89075 Ulm

Tel. Nr.: 0731 7047–0

Fax: 0731 7047–297

Figure 1

RIGHT! (because backstop opening at **back**)

WRONG! (because backstop opening at **front**)

Explanation of international symbols

-  Caution
-  Consult instructions for use
-  Keep dry
-  Temperature limitation
-  Keep away from sunlight
-  Do not use if package is damaged
-  Do not re-use
-  Do not re-sterilize
-  Does not contain any latex
-  Use by
-  **LOT** Batch number
-  **STERILE** Sterilized using steam
-   Manufacturer
-  Image of the prefilled syringe